

# **Gamma spectrometric determination of the activities of natural radionuclides**

**$\gamma$ -SPEKT/NATRAD**

Authors:

U.-K. Schkade

A. Heckel

H. Wershofen

M.-O. Aust

F. Bruchertseifer

A. Dalheimer

S. Hofmann

G. Kanisch

C. Kowalik

F. Ober

K. Rupprecht

Drafting Committee  
(Redaktionsausschuss der Messanleitungen des Bundes)

## CONTENT

<b>1</b>	<b>INTRODUCTION.....</b>	<b>1</b>
<b>2</b>	<b>SPECIAL REQUIREMENTS PLACED ON THE MEASUREMENT SYSTEM AND ITS SURROUNDINGS .....</b>	<b>1</b>
2.1	REQUIREMENTS RELATED TO THE PLACE OF INSTALLATION.....	1
2.2	REQUIREMENTS RELATED TO THE GAMMA SPECTROMETRIC MEASUREMENT SYSTEM.....	3
<b>3</b>	<b>CONSIDERATIONS ON THE RADIOACTIVE EQUILIBRIUM IN DECAY CHAINS.....</b>	<b>3</b>
<b>4</b>	<b>CONSIDERING SELECTED NATURAL RADIONUCLIDES .....</b>	<b>4</b>
4.1	URANIUM-238 .....	5
4.2	RADIUM-226 .....	7
4.3	LEAD-210.....	8
4.4	URANIUM-235 .....	9
4.5	ACTINIUM-227 .....	9
4.6	THORIUM-232.....	10
4.7	RADIUM-228 .....	11
4.8	THORIUM-228.....	11
4.9	POTASSIUM-40 .....	12
<b>ANNEX A</b>	<b>NUCLEAR DATA OF SELECTED RADIONUCLIDES .....</b>	<b>13</b>
<b>REFERENCES.....</b>		<b>19</b>

# **Gamma spectrometric determination of the activities of natural radionuclides**

## **1 Introduction**

Besides K-40 and cosmogenic radionuclides, the radionuclides of the three natural decay chains also belong to the natural radionuclides (see Figure 1). Within the scope of the Directive on Emission and Immission Monitoring regarding Mining (REI Bergbau) [1] the following nuclides

- uranium-238, radium-226 and lead-210 from the uranium-radium decay chain;
- uranium-235 and actinium-227 from the uranium-actinium decay chain, and
- thorium-232, radium-228 and thorium-228 from the thorium decay chain

are usually determined using gamma spectrometry. Some of these long-lived radionuclides can only be measured quantitatively by means of gamma spectrometry via the progeny nuclide after the radioactive equilibrium between the parent and the progeny radionuclides has been attained.

The fundamentals of gamma spectrometry are provided in the General Chapter  $\gamma$ -SPEKT/GRUNDL [2], more detailed information on background and interferences are given in the General Chapters  $\gamma$ -SPEKT/NULLEF [3] and  $\gamma$ -SPEKT/INTERF [4] of this Procedures' Manual.

## **2 Special requirements on the measurement system and its surroundings**




### **2.1 Requirements related to the place of installation**

Natural radionuclides are ubiquitous. They, therefore, occur in the construction materials used for flooring, walls and ceilings as well as in the air of laboratories and in the materials used to manufacture both the detectors and the shieldings.

In this context, particular attention is paid to the isotopes of the noble gas radon which are emanating from the above-mentioned materials into the laboratory air and contribute to the background spectrum. In particular when determining the activities of natural radionuclides, adequate measures must be taken to minimize the contents of radon isotopes and their progenies and to prevent, as far as possible, the contents to vary in order to be able to evaluate the pulse height spectrum correctly [3].

**Note:**

- High radon concentrations may occur in the basement of old buildings.
- Seasonal changes in the concentration of radon isotopes and their progenies must be taken into account, in particular when radon concentrations are low.




**Fig. 1:** The natural decay chains with decay mode and half-life as well as branching (nuclear data as of November 2024)

## 2.2 Requirements related to the gamma spectrometric measurement system

To determine very low activities of natural nuclides reliably, an ultra-low-level measurement system is usually used. Its components are made of materials whose own activity due to natural radionuclides is very low. Ideally, this measurement system should be placed in the centre of a room whose construction material exhibits activity contents of natural radionuclides as low as possible, the distance between the measurement system and the walls being as large as possible given the space of the room. The other requirements are described in [2: Section 2].

Some of the natural radionuclides can be only measured by use of gamma peaks in the low-energy region of the pulse height spectrum below 100 keV (see Figure 2). It is, therefore, an advantage if besides the current p-type germanium detectors an n-type detector can also be used.



**Fig. 2:** Pulse height spectrum of a soil sample (energy region up to 160 keV) used to determine the activity concentration of Pb-210

## 3 Considerations on the radioactive equilibrium in decay chains

Radionuclides of the natural decay chains that do not emit gamma rays can only be determined indirectly via their gamma emitting progenies. For such measurements, it

must be ensured that the radioactive equilibrium between the parent radionuclide and the progenies is reached.

In environmental samples, the radioactive equilibrium may be disturbed due to the different chemical and/or biochemical properties of the elements concerned. Examples are:

- different solubility of radium and thorium in water;
- different transfer behaviours of parent and progeny radionuclides in the *soil-vegetation-animal-milk chain*;
- disturbance of the radioactive equilibrium during the sampling (e.g. due to the use of filter materials);
- radon escaping from the sample.

The following measures are recommended:

- Storing the samples to be measured for a sufficiently long time prior to starting the measurement or keeping them for a check measurement that may be necessary. If the aim is to reach the equilibrium between a long-lived parent radionuclide and its short-lived progeny radionuclide, then the waiting time should amount to at least six half-lives of the progeny radionuclide.
- Transferring the sample material (in particular, when measuring short-lived progeny of Rn-222) into gas-tight measurement containers and waiting until the radioactive equilibrium to be reached. The dead volume between the sample filled in and the lid of the measurement cell should be as small as possible.

With respect to parent-progeny radionuclides, particular attention must be paid in order to replace the half-life of a short-lived progeny radionuclide by the half-life of the long-lived parent radionuclide, so that no physically unrealistic activities are calculated when correcting the decay to a reference point in time of an environmental sample which lies before the time of the measurement – usually, the sampling time. Then, both the radioactive decay of the respective radionuclide and growth from its parent radionuclide must be taken into account. Such examples are the radionuclide pairs: Ra-226/Pb-210, Ra-228/Th-228, U-238/Th-234.

Additional information is given in the General Chapter  $\gamma$ -SPEKT/ GRUNDL of this Procedures' Manual [2].

## 4 Considering selected natural radionuclides

In the following, the natural radionuclides required to be measured in accordance with the REI Bergbau directive will be discussed thoroughly [1]. For certain cases, the

requirement of correction for self-absorption losses and true coincidence summing will be addressed.

The tables below list the photon energies ( $E_\gamma$ ) and emission intensities ( $p_\gamma$ ) of the radio-nuclides considered as well as that of their interfering nuclides. In the case of the decay chains at equilibrium, the emission intensities stated here each refer to one decay of the parent radionuclide and already take possible branching into account. An example is the determination of Th-228 content via the progeny Tl-208 in Section 4.8. The Tl-208 emission intensities listed there already consider that the Tl-208 activity only represents 35,93 % of its parent nuclide Th-228 due to a branching of the decay chain at Bi-212 (see Figure 1).

**Note:**

Nuclear data of the regarded radionuclides refer to the sources mentioned in the General Chapter KERNDATEN of this Procedures' Manual [5].

## 4.1 Uranium-238


With its half-life of  $4,468 \cdot 10^9$  years, uranium-238 (U-238) is the parent radionuclide in the uranium-radium decay chain. Using gamma spectrometry, U-238 can only be determined via its progenies Th-234 and Pa-234m.

In most cases, the 63,3 keV peak of Th-234 is used for this purpose. At such a low energy, the matrix differences between the calibration and counting source significantly influence the results due to specific self-absorption effects and must be considered correspondingly [2]. In addition, the peak at 63,3 keV is subject to an interference caused by the peak of Th-232 at 63,8 keV with an emission intensity of 0,259 %.

**Note:**

To determine the U-238 content, the doublet peak of Th-234 at the energies of 92,4 keV and 92,8 keV with the combined emission intensity of 4,33 % (2,18 % + 2,15 %), should not be used.

Due their higher energies, the peaks of Pa-234m at 1001,0 keV and 766,4 keV are better suited to evaluate the pulse height spectrum. Their low emission intensities, however, lead to poor detection limits (see Figure 3).



**Fig. 3:** Evaluating the Pa-234m peak in the pulse height spectrum of a soil sample used to determine of the specific activity of U-238

The peaks suitable for the evaluation are compiled in Table 1.

**Tab. 1:** Nuclear data of the radionuclides to be determined and their interfering nuclides when determining the U-238 content by means of gamma spectrometry

Radionuclide to be determined	$E_\gamma$	$p_\gamma$	Interfering radionuclide	$E_\gamma$	$p_\gamma$
	in keV	in %		in keV	in %
Th-234	63,30	3,75	Th-232	63,81	0,259
	92,38	2,18	Th-234 ( $K\alpha_1 + K\alpha_2$ )	93,35	5,6
	92,80	2,15			
Pa-234m	766,36	0,323			
	1001,03	0,847			

In the case of samples of natural origin, the nearly constant natural U-238/U-235 activity ratio of 21,7 may be referred to for plausibility reasons (see Section 4.4).

## 4.2 Radium-226

With a half-life of 1600 years, radium-226 (Ra-226) is a long-lived radionuclide in the uranium-radium decay chain. To determine this radionuclide, two approaches are possible:

- evaluation of the 186,2 keV peak of Ra-226 in the pulse height spectrum, or
- evaluation of the gamma peaks of the short-lived progeny radionuclides of Rn-222, Pb-214 and Bi-214, in the pulse height spectrum after reaching the radioactive equilibrium between Ra-226, Rn-222, Pb-214 and Bi-214 (see Figure 4).


If the Ra-226 peak at 186,2 keV is evaluated, its overlapping with the U-235 peak at 185,7 keV must be considered. In natural sample materials, the activity content of U-235 amounts to approx.  $1/22$  of the 1 of Ra-226. In comparison, the emission intensity of the U-235 peak at 185,7 keV is by a factor of 16 higher than the emission intensity of the Ra-226 line at 186,2 keV.

Assuming that the sample is at radioactive equilibrium, this overlapping is not problematic. This means that the contents of U-238 and Ra-226 are the same and therefore, the content of U-235 is also known. Environmental samples, however, are not always in radioactive equilibrium, while industrial residues and products containing naturally occurring radioactive material (NORM) are never in equilibrium. Therefore, the correction due to the overlapping of the Ra-226 peak at 186,2 keV with the interfering U-235 peak at 185,7 keV is only possible if the U-235 activity is determined based on other gamma peaks or if it is derived from the U-238 activity.

**Note:**

The U-238 content of a sample can also be obtained by means of other methods such as alpha spectrometry, X-ray fluorescence measurement or mass spectrometry.

Gas-tight measurement containers are required to determine the activities of the short-lived progeny radionuclides of Rn-222, since otherwise, the pulse counts acquired for Pb-214 and Bi-214 in the pulse height spectrum are too low due to Rn-222 emanating from the measurement container. The period of time between the filling of the sample into the gas-tight measurement container and the beginning of the measurement should amount to at least 23 days, this period being determined by the half-life of Rn-222 which amounts to 3,82 days. A summation correction is necessary when evaluating the pulse height spectrum [2, 6].



**Fig. 4:** Evaluating the Pb-214 and the Bi-214 peaks in the pulse height spectrum of a soil sample used to determine the specific activity of Ra-226

A few peaks suitable for the evaluation are compiled in Table 2.

**Tab. 2:** Nuclear data of the radionuclides to be determined and their interfering nuclides when determining the Ra-226 content by means of gamma spectrometry

Radionuclide to be determined	$E_\gamma$	$p_\gamma$	Interfering radionuclide	$E_\gamma$	$p_\gamma$
	in keV	in %		in keV	in %
Ra-226	186,21	3,56	U-235	185,72	57,0
Pb-214	295,22	18,41			
	351,93	35,6	Bi-211	351,03	13,00
Bi-214	609,31	45,49			
	1120,29	14,91			
	1764,49	15,31			

### 4.3 Lead-210

Lead-210 (Pb-210) is a radionuclide in the uranium-radium decay chain; its half-life is 22,23 years. The Pb-210 content is determined directly via the gamma peak at the relatively low energy of 46,5 keV (see Table 3). At this energy, the differences in the material composition of the calibration and counting source play a key role due to self-absorption; adequate correction factors must be applied [2].

**Tab. 3:** Nuclear data of the radionuclide to be evaluated for the gammaspectrometric determination of the Pb-210 content

Radionuclide to be determined	$E_\gamma$	$p_\gamma$	Interfering radionuclide	$E_\gamma$	$p_\gamma$
	in keV	in %		in keV	in %
Pb-210	46,54	4,25			

#### 4.4 Uranium-235

Uranium-235 (U-235) is the parent radionuclide of the uranium-actinium decay chain; its half-life is  $7,04 \cdot 10^8$  years. The U-235 content can be determined directly by means of the gamma peaks listed in Table 4 at 143,8 keV, 163,4 keV, 185,7 keV and 205,3 keV. The following aspects must be taken into account when evaluating these gamma peaks in the pulse height spectrum:


- To evaluate the gamma peak at 185,7 keV, which has the highest emission intensity, the activity of Ra-226 must be known (see Section 4.2).
- The Ra-223 gamma peak at 144,2 keV interferes with the gamma peak at 143,8 keV (see Figure 2).
- The two other gamma peaks have emission intensities that are so low that they are usually not considered to evaluate samples with low U-235 content.

**Tab. 4:** Nuclear data of the radionuclide considered and its interfering nuclides when determining the U-235 content by means of gamma spectrometry

Radionuclide to be determined	$E_\gamma$	$p_\gamma$	Interfering radionuclide	$E_\gamma$	$p_\gamma$
	in keV	in %		in keV	in %
U-235	143,77	10,94	Ra-223	144,27	3,36
	163,36	5,08			
	185,72	57,0	Ra-226	186,21	3,56
	205,32	5,02			

#### 4.5 Actinium-227

Actinium-227 (Ac-227) is a radionuclide in the uranium-actinium decay chain; its half-life is 21,77 years. Ac-227 can only be determined by gamma spectrometry using its short-lived progenies Th-227, Ra-223 and Rn-219. Preferably, the Th-227 peak at 236,0 keV in the pulse height spectrum is evaluated (see Section 5 and Table 5). The other potentially evaluable peaks either have lower emission intensities or the gamma peaks of other radionuclides interfere with them.



**Fig. 5:** Evaluation of the 236,0 keV peak of Th-227 in the pulse height spectrum of a sponge sample (Greenland) used to determine the specific activity of Ac-227

**Tab. 5:** Nuclear data of the radionuclide to be evaluated when determining the Ac-227 content by means of gamma spectrometry

Radionuclide to be determined	$E_\gamma$	$p_\gamma$	Interfering radionuclide	$E_\gamma$	$p_\gamma$
	in keV	in %		in keV	in %
Th-227	235,96	12,6			

## 4.6 Thorium-232

With a half-life of  $1,40 \cdot 10^{10}$  years, thorium-232 (Th-232) is the parent radionuclide of the thorium decay chain. At 63,8 keV, Th-232 has a gamma peak with a very low emission intensity of 0,259 %. A gamma peak of Th-234 at 63,3 keV with a higher emission intensity of 3,75 % interferes with this gamma peak (see Table 1), so that Th-232 cannot be determined directly by means of gamma spectrometry in environmental samples.

Determining Th-232 via its progenies Ac-228, Pb-212 and Tl-208 is only possible if these radionuclides are at radioactive equilibrium with each other and with their parent radionuclide Th-232. Environmental samples, however, are not always in radioactive equilibrium, while industrial residues and products containing naturally occurring radioactive material (NORM) are never in equilibrium. This is caused by the radionuclide Ra-228, which is located between Th-232 and Ac-228 in the decay chain, and can modify the radioactive equilibrium due to its chemical properties and mobility.

## 4.7 Radium-228

Radium-228 (Ra-228) is a radionuclide of the thorium decay chain; its half-life is 5,75 years. The Ra-228 content is determined by evaluating the gamma peaks of its progeny Ac-228, which is at radioactive equilibrium with Ra-228 after approx. 30 hours due to its short half-life of 6,15 hours. Here, a coincidence summation correction must always be performed [2]. In Table 6, gamma peaks are listed which are suitable for the evaluation.

If the activities of Ra-224, Pb-212 and Tl-208 (whose gamma peaks are listed in Table 7) do not significantly deviate from the activity of Ac-228, then these radionuclides may also be used to determine the activity of Ra-228.

**Tab. 6:** Nuclear data of the radionuclide to be evaluated and of the interfering nuclide when determining the Ra-228 content by means of gamma spectrometry

Radionuclide to be determined	$E_\gamma$	$p_\gamma$	Interfering radionuclide	$E_\gamma$	$p_\gamma$
	in keV	in %		in keV	in %
Ac-228	209,25	3,97	Ra-223	338,28	2,85
	338,32	11,40			
	911,20	26,20			
	968,96	15,90			

## 4.8 Thorium-228

Thorium-228 (Th-228) is a radionuclide of the thorium decay chain; its half-life is 1,91 years. At radioactive equilibrium, Th-228 content may be determined by evaluating the gamma peaks of its short-lived progenies Ra-224, Pb-212 and Tl-208. Regarding radioactive equilibrium, it is important to differentiate between solid and liquid samples:

- In liquid samples, radioactive equilibrium is certainly reached after a period of 23 days.
- In solid samples, radioactive equilibrium is prevalent at the time of sampling.

When evaluating the gamma peak of Ra-224 at 241,0 keV, the interfering gamma peak of Pb-214 at an energy of 242,0 keV must be taken into account (see Figure 5).

If the content of Th-228 is determined via the gamma peaks of the radionuclides Pb-212 and Tl-208, Rn-220 must remain in the sample matrix which is usually given due to its low half-life of 55,8 seconds. Preferably, the gamma peak of Pb-212 at 238,6 keV is used for the evaluation where coincidence summation corrections must always be applied [2].

**Tab. 7:** Nuclear data of the radionuclides to be evaluated and the interfering nuclides when determining the Th-228 content by means of gamma spectrometry

Radionuclide to be determined	$E_\gamma$	$p_\gamma$	Interfering radionuclide	$E_\gamma$	$p_\gamma$
	in keV	in %		in keV	in %
Ra-224	240,99	4,12	Pb-214	242,00	7,27
Pb-212	238,63	43,60			
	300,09	3,18	Th-227	300,50	0,014
			Pa-231	300,06	2,41
TI-208	277,37	2,37	Ac-228	278,80	0,235
	583,19	30,54	Ac-228	583,39	0,120
	860,53	4,46			
	2614,51	35,84			

#### 4.9 Potassium-40

Potassium-40 (K-40) is present in the isotopic mixture of natural potassium (K-39, K-40, K-41) at the rate of 0,0117 weight by weight. It decays with a half-life of  $1,25 \cdot 10^9$  years to Ca-40 ( $\beta^-$  decay) and Ar-40 (electron capture), respectively. K-40 is determined via its 1460,8 keV gamma peak (see Figure 4).

The interference due to Ac-228 shown in Table 8 must be taken into account for example when determining activity content in sewage sludge and/or in various NORM substances in which – contrary to environmental samples – elevated activity contents of Ra-228 occur as compared to K-40.

**Tab. 8:** Nuclear data of the radionuclide to be evaluated and of the interfering nuclide when determining the K-40 content by means of gamma spectrometry

Radionuclide to be determined	$E_\gamma$	$p_\gamma$	Interfering radionuclide	$E_\gamma$	$p_\gamma$
	in keV	in %		in keV	in %
K-40	1460,82	10,55	Ac-228	1459,13	0,87

## Annex A

### Nuclear data of selected radionuclides

Sources of nuclear decay data are compiled in the General Chapter KERNDATEN of this Procedures' Manual [5]. Table A1 compiles the gamma photon energies in the energy region above 25 keV and the associated emission intensities corrected for branching for selected natural radionuclides. The listed emission intensities of X-ray peaks are not corrected for branching, which is indicated by an asterisk (\*) following the emission intensity.

**Note:**

The emission intensities of the characteristic X-ray photons apply to the corresponding nuclear transformations. It must be pointed out that characteristic X-rays may also occur in the form of fluorescence radiation. When determining the activity of natural radionuclides in environmental samples, the contribution of fluorescence radiation to the net count rates is, however, negligible in most instances.

In addition, the column "Other" contains artificial radionuclides, which are generally used in calibration standards, and activation products, which may be generated in the measurement system by neutrons originating in cosmic radiation [3: Section 2, 5].

**Tab. A1:** Photon energies and emission intensities of selected radionuclides

Peak $E_\gamma$ in keV	Th-232 and progeny radionuclides		Ra-226 and progeny radionuclides		U-238/U-235 and progeny radionuclides		Ac-227 and progeny radionuclides		Other	
	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$ or origin
25,64					Th-231	0,139				
39,86	Bi-212	0,00107								
46,54			Pb-210	0,0425						
50,13										
53,16										
53,20					U-234	0,00125				
53,23			Pb-214	0,0106						
53,47										
57,75	Ac-228	0,0047								
59,54										
63,30					Th-234	0,0375				
63,81	Th-232	0,00259								
66,75										
67,67					Th-230	0,00377				
72,80	Pb-K $\alpha$ 2	0,0077*								
72,87										
74,82	Bi-K $\alpha$ 2	0,107*	Bi-K $\alpha$ 2	0,0626*					Tl-K $\alpha$ 1	0,01225*
74,97	Pb-K $\alpha$ 1	0,0361*								

Peak in keV	Th-232 and progeny radionuclides		Ra-226 and progeny radionuclides		U-238/U-235 and progeny radionuclides		Ac-227 and progeny radionuclides		Other		
	$E_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$ or origin
76,86				Po-K $\alpha$ 2	0,00426*						
77,11	Bi-K $\alpha$ 1	0,177*		Po-K $\alpha$ 1	0,00710*						
79,29										Ba-133	0,0263
79,61							Th-227	0,0189			
79,72										Ba-133	0,333
81,00							Rn-K $\alpha$ 2	0,150*			
81,07	Rn-K $\alpha$ 2	0,0013	Rn-K $\alpha$ 2	0,00192*			Rn-K $\alpha$ 1	0,249*			
83,78	Rn-K $\alpha$ 1	0,00215	Rn-K $\alpha$ 1	0,00317*							
84,21					Th-231	0,0670	Ra-K $\alpha$ 2	0,0017*			
84,37	Th-228	0,0119					Ra-K $\alpha$ 1	0,0028*		Cd-109	0,0366
84,94	Pb-K $\alpha$ 1	0,0117*		Po-K $\alpha$ 1	0,00245*						
85,43						Th-234	0,0218				
87,34	Bi-K $\alpha$ 1	0,0612*				Th-234	0,0215				
88,03							Th-227	0,0137			
88,47							Rn-K $\alpha$ 1	0,0870*			
89,81										Se-75	0,0335
89,95	Th-K $\alpha$ 2	0,025*									
92,38											
92,80											
93,35	Th-K $\alpha$ 1	0,041*									
93,93											
94,86											
96,73											
97,90											
99,51	Ac-228	0,0126									
102,27						Th-231	0,00441				
105,31											
105,60	Th-K $\alpha$ 1	0,00155*									
109,19											
109,70											
112,81											
115,18	Pb-212	0,00624									
121,12											
122,06											
122,32											
129,07	Ac-228	0,0250									
136,00											
136,47											
139,20											
143,77											
144,27											
145,44											
145,84	Ac-228	0,00169									
153,97	Ac-228	0,00754									

Peak $E_\gamma$ in keV	Th-232 and progeny radionuclides		Ra-226 and progeny radionuclides		U-238/U-235 and progeny radionuclides		Ac-227 and progeny radionuclides		Other	
	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$ or origin
154,21							Ra-223	0,0584		
159,50							Ge-77m	(n, $\gamma$ )		
160,61							Ba-133	0,00638		
163,36					U-235	0,0508				
174,88					U-235	0,570			Ge-71m	(n, $\gamma$ )
185,72									Cu-66	(n, $\gamma$ )
186,01			Ra-226	0,0356					Ge-70m	(n, $\gamma$ )
186,20										
191,35	Ac-228	0,00133							Ge-71m	(n, $\gamma$ )
197,90									Se-75	0,0146
198,61										
202,12					U-235	0,0108				
205,32					U-235	0,0502				
209,25	Ac-228	0,0397					Th-227	0,0111		
210,65									Ge-77	(n, $\gamma$ )
215,51									Ba-133	0,00450
223,24							Th-227	0,126		
235,96										
238,63	Pb-212	0,436					Th-227	0,070		
240,99	Ra-224	0,0412								
242,00			Pb-214	0,0727						
256,23										
258,87			Pb-214	0,00532						
264,66									Se-75	0,588
269,46										
270,25	Ac-228	0,0355								
271,23										
274,80			Pb-214	0,00362						
276,40										
277,37	Tl-208	0,0237								
278,24										
278,80	Ac-228	0,00235								
279,20										
279,54										
283,69					Pa-231	0,0165				
286,12							Th-227	0,0154		
288,18	Bi-212	0,0032								
295,22			Pb-214	0,184						
300,05										
300,06										
300,09	Pb-212	0,0318								
302,67										
302,85										
303,92										

Peak $E_\gamma$ in keV	Th-232 and progeny radionuclides		Ra-226 and progeny radionuclides		U-238/U-235 and progeny radionuclides		Ac-227 and progeny radionuclides		Other	
	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$ or origin
304,52							Th-227	0,012		
320,08									Cr-51	0,0989
321,65	Ac-228	0,00232					Ra-223	0,0406		
323,87							Th-227	0,027		
328,04	Bi-212	0,00121					Th-227	0,0105		
328,00	Ac-228	0,0304					Ra-223	0,0285		
329,85										
330,04					Pa-231	0,0136	Bi-211	0,130		
332,37	Ac-228	0,0037							Ba-133	0,621
334,38									Ba-133	0,089
338,28									Se-75	0,114
338,32	Ac-228	0,1140								
340,97	Ac-228	0,00405					Rn-219	0,0675		
351,03							Pb-211	0,0383		
351,93			Pb-214	0,356			Pb-211	0,0181		
356,01							Ra-223	0,0128		
383,85									Cs-134	0,0148
386,77			Bi-214	0,00296					Be-7	0,104
388,88			Bi-214	0,00394						
400,66										
401,81										
404,83										
409,46	Ac-228	0,0202								
427,15										
445,03										
452,98	Bi-212	0,00340								
454,77			Bi-214	0,00288						
463,00	Ac-228	0,0445								
475,37										
477,60										
480,43			Pb-214	0,00337						
487,09			Pb-214	0,00433						
503,82	Ac-228	0,00171								
508,96	Ac-228	0,00510								
510,74	Tl-208	0,0808								
558,46									Cd-114	(n, $\gamma$ )
562,50	Ac-228	0,00890							Cs-134	0,0834
563,25									Cs-134	0,154
569,33										
570,88	Ac-228	0,00190								
580,13			Pb-214	0,00369						
583,19	Tl-208	0,305								
583,39	Ac-228	0,00120								
595,85									Ge-74	(n,n', $\gamma$ )

Peak $E_\gamma$ in keV	Th-232 and progeny radionuclides		Ra-226 and progeny radionuclides		U-238/U-235 and progeny radionuclides		Ac-227 and progeny radionuclides		Other	
	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$ or origin
600,70									Ge-74	(n,n', $\gamma$ )
604,72									Cs-134	0,976
608,35									Ge-74	(n,n', $\gamma$ )
609,31			Bi-214	0,455					Cd-114	(n, $\gamma$ )
651,26									Cs-137	0,850
661,66									Cu-63	(n,n', $\gamma$ )
665,45			Bi-214	0,0153					Ge-72	(n,n', $\gamma$ )
669,60										
693,40										
703,11			Bi-214	0,00479						
719,86			Bi-214	0,00393						
726,88	Ac-228	0,00680								
727,33	Bi-212	0,0665								
755,31	Ac-228	0,0103								
763,45	Tl-208	0,0647								
765,80									Nb-95	0,998
766,36										
768,36			Bi-214	0,0489						
772,29	Ac-228	0,0152								
782,14	Ac-228	0,00500								
785,37	Bi-212	0,0111								
785,96			Pb-214	0,0106						
794,94	Ac-228	0,0431								
795,86									Cs-134	0,855
801,95									Cs-134	0,0869
803,00									Pb-206	(n,n', $\gamma$ )
805,89									Cd-114	(n, $\gamma$ )
806,17			Bi-214	0,0126						
830,48	Ac-228	0,00610								
831,98										
834,85									Mn-54	0,9998
835,70	Ac-228	0,0170								
839,04			Pb-214	0,00587						
840,37	Ac-228	0,00970								
860,53	Tl-208	0,0446								
867,90										
893,41	Bi-212	0,0038								
904,20	Ac-228	0,0078								
911,20	Ac-228	0,262								
934,06			Bi-214	0,0310						
962,00										
964,79	Ac-228	0,0499								
968,96	Ac-228	0,159								
1001,03										
					Pa-234m	0,00847				

Peak $E_\gamma$ in keV	Th-232 and progeny radionuclides		Ra-226 and progeny radionuclides		U-238/U-235 and progeny radionuclides		Ac-227 and progeny radionuclides		Other	
	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$ or origin
1038,61									Cs-134	0,00991
1051,96			Bi-214	0,00324					Pb-207	(n,n', $\gamma$ )
1063,60									Cu-65	(n,n', $\gamma$ )
1078,63	Bi-212	0,00550	Bi-214	0,149					Ge-(75)	(n, $\gamma$ )
1115,54			Bi-214	0,0164					Cs-134	0,0179
1120,29			Bi-214	0,00454					Co-60	0,9985
1127,80			Bi-214	0,0583					Cu-63	(n,n', $\gamma$ )
1155,19			Bi-214	0,0144					Co-60	0,9998
1167,97			Bi-214	0,0397					Cs-134	0,0302
1173,23			Bi-214	0,00795					Cu-63	(n,n', $\gamma$ )
1207,68			Bi-214	0,0133					K-40	0,106
1238,11			Bi-214	0,0239						
1247,10	Ac-228	0,00524	Bi-214	0,00707						
1280,96			Bi-214	0,00322						
1327,00			Bi-214	0,00401						
1332,49			Bi-214	0,0302						
1365,19			Bi-214	0,00707						
1377,67			Bi-214	0,00322						
1385,31			Bi-214	0,0116						
1401,50			Bi-214	0,0491						
1407,98			Bi-214	0,0151						
1412,00			Bi-214	0,0152						
1459,13	Ac-228	0,00870	Bi-214	0,0152						
1460,82			Bi-214	0,0290						
1495,90	Ac-228	0,00920	Bi-214	0,0213						
1501,59	Ac-228	0,00513	Bi-214	0,0401						
1509,23			Bi-214	0,0302						
1512,70	Bi-212	0,00290	Bi-214	0,00707						
1538,50			Bi-214	0,00707						
1543,32			Bi-214	0,0116						
1580,53	Ac-228	0,00620	Bi-214	0,0491						
1583,22			Bi-214	0,0151						
1588,20	Ac-228	0,0306	Bi-214	0,0152						
1599,31			Bi-214	0,0290						
1620,74	Bi-212	0,0151	Bi-214	0,0284						
1630,62	Ac-228	0,0152	Bi-214	0,1531						
1638,27	Ac-228	0,0046	Bi-214	0,0034						
1661,28			Bi-214	0,0203						
1729,60			Bi-214	0,0116						
1764,49			Bi-214	0,0406						
1838,36			Bi-214	0,0116						
1847,42			Bi-214	0,0491						
2118,55			Bi-214	0,0151						
2204,21			Bi-214	0,0290						

Peak $E_\gamma$ in keV	Th-232 and progeny radionuclides		Ra-226 and progeny radionuclides		U-238/U-235 and progeny radionuclides		Ac-227 and progeny radionuclides		Other	
	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$	Nuclide	$p_\gamma$ or origin
2293,40			Bi-214	0,0031						
2447,86			Bi-214	0,0155						
2614,51	Tl-208	0,358								

## References

- [1] Richtlinie zur Emissions- und Immissionsüberwachung bei bergbaulichen Tätigkeiten (REI-Bergbau) from 10.11.2020. Available at: <https://www.bfs.de/SharedDocs/Downloads/BfS/DE/fachinfo/ion-altlasten/rei-bergbau.pdf?blob=publicationFile&v=3>. [Last access 21.08.2024].
- [2] Arnold, D., Debertin, K., Heckel, A., Kanisch, G., Wershofen, H., Wilhem, C.: *Fundamentals of gamma spectrometry*.  $\gamma$ -SPEKT/GRUNDL, Version March 2018. In: Federal Ministry for the Environment, Nature Conservation, Nuclear Safety and Consumer Protection (ed.): Procedures' Manual for monitoring of radioactive substances in the environment and of external radiation. ISSN 1865-8725. Available at: <https://www.bmuv.de/WS1517-1>. [Last access 27.10.2023].
- [3] Arnold, D., Heckel, A., Wershofen, H.: *Background contributions in gamma spectrometry*.  $\gamma$ -SPEKT/NULLEF, Version June 2019. In: Federal Ministry for the Environment, Nature Conservation, Nuclear Safety and Consumer Protection (ed.): Procedures' Manual for monitoring of radioactive substances in the environment and of external radiation. ISSN 1865-8725. Available at: <https://www.bmuv.de/WS1517-1>. [Last access 21.08.2024].
- [4] Kanisch, G., Mundschenk, H., Schkade, U.-K., Wershofen, H.: *Interferences in gamma spectrometry*.  $\gamma$ -SPEKT/INTERF, Version March 2019. In: Federal Ministry for the Environment, Nature Conservation, Nuclear Safety and Consumer Protection (ed.): Procedures' Manual for monitoring of radioactive substances in the environment and of external radiation. ISSN 1865-8725. Available at: <https://www.bmuv.de/WS1517-1>. [Last access 21.08.2024].
- [5] Heckel, A., Wershofen, H.: *Nuclear Data*. KERNDATEN, Version November 2024. In: Federal Ministry for Environment, Nature Conservation, Nuclear Safety and Consumer Protection, (ed.): Procedures' Manual for moitoring of radioactive substances in the environment and of external radiation. ISSN 1865-8725. Available at: <https://www.bmuv.de/WS1517-1>. [Last access 22.08.2024].
- [6] Arnold, D., Kanisch, G., Wershofen, H.: *Summation peaks and escape peaks of selected radionuclides in gamma spectrometry*.  $\gamma$ -SPEKT/SUMESC, Version November 2018. In: Federal Ministry for the Environment, Nature Conservation, Nuclear Safety and Consumer Protection (ed.): Procedures' Manual for monitoring of radioactive

substances in the environment and of external radiation. ISSN 1865-8725. Available at: <https://www.bmuv.de/WS1517-1>. [Last access 21.08.2024].